Bruksanvisning

extendSURE™ HbA1c vätskekontroller

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Referenser

1. IFCC standardization of HbA1c. <https://www.ngsp.org/ifccngsp.asp>
2. Jeppsson J-O, Kobold U, Barr J, Finke A, Hoebel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med 2002; 40(1): 78-89.
3. Hoelzel W, Weykamp C, Jeppsson J-O, Miedema, Barr JR, Goodall I, Hoshino T, John, WG, Kobold U, Little R, Mosca A, Mauri P, Paroni R, Susanto F, Takei I, Theinpoint L, Umemoto M, Wiedmeyer H-M. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardisation Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50(1): 166-174.
4. Weykamp, C. et al. The IFCC Reference Measurement System for HbA1c: A 6-Year Progress Report Clinical Chemistry 54:2: 240-248 (2008).
5. Andrea Geistanger, Sabine Arends, Christoph Berding, Tadao Hoshino, Jan-Olof Jeppsson, Randie Little, Carla Siebelder and Cas Weykamp on behalf of the IFCC Working Group on Standardization of HbA1c: Statistical Methods for Monitoring the Relationship between the IFCC Reference Measurement Procedure for Hemoglobin A1c and the Designated Comparison Methods in the United States, Japan and Sweden. Clin Chem 2008, 54 (8): 1379-85.
6. Andrea Mosca, Ian Goodall, Tadao Hoshino, Jan O. Jeppsson, W. Garry John, Randie R. Little, Kor Miedema, Gary L. Myers, Hans Reinauer, David B. Sacks and Cas W. Weykamp. Global standardization of glycosylated hemoglobin measurement: the position of the IFCC Working Group. Clin Chem Lab Med 2007, 45(8): 1077-1080.
7. Little RR, Rohlfing CL, Sacks DB. Status of HbA1c measurement and goals for improvement: From chaos to order for improving diabetes care. Clin Chem 2011; 57:205-214.
8. National Glycohemoglobin Standardization Program. NGSP protocol. (Uppdaterat januari 2019) <http://www.ngsp.org/docs/Protocol.pdf> (Hämtat november 2021).

71 Whiteleigh Avenue, Christchurch 8011, NYA ZEELAND Tel. +64 3 343 3345

EC REP Emergo Europe
Westervoortsedijk 60,
6827 AT Arnhem,
Nederländerna

CH REP Beckman Coulter International S.A.
Rue Juste-Olivier 22,
1260 – Nyon 1,
Schweiz