Upute za uporabu

extendSURE™ - tekuće kontrole HbA1c

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

Reference:

1. IFCC standardization of HbA1c (IFCC standardizacija HbA1c). <https://www.ngsp.org/ifccngsp.asp>
2. Jeppsson J-O, Kobold U, Barr J, Finke A, Hoebel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the Measurement of HbA1c in Human Blood (Odobrena IFCC referentna metoda za mjerenje HbA1c u ljudskoj krvi). Clin Chem Lab Med 2002; 40(1): 78-89.
3. Hoelzel W, Weykamp C, Jeppsson J-O, Miedema, Barr JR, Goodall I, Hoshino T, John, WG, Kobold U, Little R, Mosca A, Mauri P, Paroni R, Susanto F, Takei I, Theinpoint L, Umemoto M, Wiedmeyer H-M. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardisation Schemes in the United States, Japan, and Sweden (IFCC referentni sustav za mjerenje hemoglobina A1c u ljudskoj krvi i nacionalni programi standardiziranja u Sjedinjenim Državama, Japanu i Švedskoj): A Method-Comparison Study (studija uspoređivanja metoda). Clin Chem 2004; 50(1): 166-174.
4. Weykamp, C. et al. The IFCC Reference Measurement System for HbA1c (IFCC referentni mjerni sustav za HbA1c): A 6-Year Progress Report Clinical Chemistry 54:2 (6-godišnje izvješće o napretku - klinička kemija 54:2): 240-248 (2008).
5. Andrea Geistanger, Sabine Arends, Christoph Berding, Tadao Hoshino, Jan-Olof Jeppsson, Randie Little, Carla Siebelder i Cas Weykamp u ime IFCC radne skupine za standardiziranje HbA1c: Statistical Methods for Monitoring the Relationship between the IFCC Reference Measurement Procedure for Hemoglobin A1c and the Designated Comparison Methods in the United States, Japan and Sweden (Statističke metode za praćenje veze između IFCC referentne mjerne procedure hemoglobina A1c i metode za usporedbu u Sjedinjenim Državama, Japanu i Švedskoj). Clin Chem 2008, 54; 8(1): 1379-85.
6. Andrea Mosca, Ian Goodall, Tadao Hoshino, Jan O. Jeppsson, W. Garry John, Randie R. Little, Kor Miedema, Gary L. Myers, Hans Reinauer, David B. Sacks i Cas W. Weykamp. Global standardization of glycated hemoglobin measurement: the position of the IFCC Working Group (Globalno standardiziranje mjerenja glikiranog hemoglobina: stajalište IFCC radne skupine). Clin Chem Lab Med 2007, 45; 8(1): 1077-1080.
7. Little RR, Rohlfing CL, Sacks DB. Status of HbA1c measurement and goals for improvement (Status mjerenja HbA1c i ciljevi koji se tiču poboljšanja): From chaos to order for improving diabetes care (Od kaosa do reda u poboljšanju skrbi o oboljelima od dijabetesa). Clin Chem 2011; 57:205-214.
8. National Glycohemoglobin Standardization Program (Nacionalni program standardiziranja glikohemoglobina). NGSP protocol (NGSP protokol). (Ažuriran u siječnju 2019.) <http://www.ngsp.org/docs/Protocol.pdf> (Pristup u studenom 2021.).

71 Whiteleigh Avenue, Christchurch 8011, NOVI ZELAND Telefon +64 3 343 3345

EC REP

Emergo Europe
Westervoortsedijk 60,
6827 AT Arnhem,
Nizozemska

CH REP

Beckman Coulter International S.A.
Rue Juste-Olivier 22,
1260 – Nyon 1,
Švicarska