

# extendSURE™

CONTROL

## LOT SPECIFIC VALUE ASSIGNMENT CARD

## LEVEL 1

|                   | UNITS    | MEAN  | RANGE       |
|-------------------|----------|-------|-------------|
| Haemoglobin A1c   | %        | 6.0   | 5.0-7.0     |
|                   | mmol/mol | 42    | 35-49       |
|                   | g/dL     | 0.57  | 0.47-0.67   |
|                   | mmol/L   | 0.35  | 0.29-0.41   |
| Total Haemoglobin | g/dL     | 13.56 | 11.18-15.94 |
|                   | mmol/L   | 8.41  | 6.94-9.88   |

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

**BECKMAN COULTER DxC CHEMISTRY SYSTEMS**  
For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT 4266



EC REP

Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New Zealand



USE BY 2024-12-31

Version 10. Revision date 07/19



IVD

EMERGO EUROPE  
Prinsesgracht 20  
2514 AP The Hague  
The Netherlands

# extendSURE™

CONTROL

## LOT SPECIFIC VALUE ASSIGNMENT CARD

## LEVEL 2

|                   | UNITS    | MEAN  | RANGE       |
|-------------------|----------|-------|-------------|
| Haemoglobin A1c   | %        | 10.7  | 8.9-12.5    |
|                   | mmol/mol | 93    | 77-109      |
|                   | g/dL     | 1.26  | 1.04-1.48   |
|                   | mmol/L   | 0.78  | 0.64-0.92   |
| Total Haemoglobin | g/dL     | 13.45 | 11.03-15.87 |
|                   | mmol/L   | 8.34  | 6.84-9.84   |

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

**BECKMAN COULTER DxC CHEMISTRY SYSTEMS**  
For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT 4266



EC REP

Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New Zealand



USE BY 2024-12-31

Version 10. Revision date 10/19



IVD

EMERGO EUROPE  
Prinsesgracht 20  
2514 AP The Hague  
The Netherlands

Printed to either side of  
225g card