

extendSURE™ CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 1

Haemoglobin A1c	UNITS	MEAN	RANGE
NGSP ALIGNED	%	5.8	4.8-6.8
IFCC ALIGNED	mmol/mol	40	33-47

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT 4266

USE BY 2024-12-31

Ver 02. Revision date 07/19



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

extendSURE™ CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 2

Haemoglobin A1c	UNITS	MEAN	RANGE
NGSP ALIGNED	%	10.2	8.5-11.9
IFCC ALIGNED	mmol/mol	88	74-102

Printed to either side of
225g card

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT 4266

USE BY 2024-12-31

Ver 02. Revision date 07/19



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands