

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 1

Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

5.7

(4.7-6.7)

IFCC ALIGNED

mmol/mol

38

(31-45)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT 4204



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

USE BY 2022-03-31



EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Ver 02. Revision date 07/19

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 2

Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

9.9

(8.2-11.6)

IFCC ALIGNED

mmol/mol

85

(71-99)

Printed to either side of
225g card

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT 4204



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

USE BY 2022-03-31



EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Ver 02. Revision date 07/19