

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 1

Haemoglobin A1c

NGSP ALIGNED

IFCC ALIGNED

UNITS

%

mmol/mol

g/dL

mmol/L

MEAN

5.8

40

0.56

0.35

13.82

mmol/L

RANGE

(4.8-6.8)

(33-47)

(0.46-0.66)

(0.29-0.41)

(11.44-16.20)

(7.10-10.04)

Total Haemoglobin

g/dL

13.82

(11.44-16.20)

mmol/L

8.57

(7.10-10.04)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4307



USE BY 2026-04-30

Version 12. Revision date 11/23



IVD



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EC REP

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 2

Haemoglobin A1c

NGSP ALIGNED

IFCC ALIGNED

UNITS

%

mmol/mol

g/dL

mmol/L

MEAN

10.6

92

1.28

0.79

RANGE

(8.8-12.4)

(76-108)

(1.06-1.50)

(0.65-0.93)

Total Haemoglobin

g/dL

mmol/L

13.86

8.60

(11.44-16.28)

(7.10-10.10)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER **DxC** CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4307



USE BY 2026-04-30

Version 12. Revision date 11/23



IVD



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EC REP

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands