

# extendSURE™

## CONTROL

### LOT SPECIFIC VALUE ASSIGNMENT CARD

### LEVEL 1

## Haemoglobin A1c

NGSP ALIGNED

IFCC ALIGNED

UNITS

%

mmol/mol

g/dL

mmol/L

MEAN

5.8

40

0.58

0.36

RANGE

(4.8-6.8)

(34-46)

(0.48-0.68)

(0.30-0.42)

## Total Haemoglobin

g/dL

mmol/L

14.58

9.04

(11.97-17.19)

(7.42-10.66)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

## BECKMAN COULTER AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Reagent Kit Part Number B00389 only

LOT

4307

CE 2797

EC REP



USE BY 2026-04-30



IVD

CH REP



Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New Zealand

EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands

Beckman Coulter International S.A.  
Rue Juste-Olivier 22,  
1260 - Nyon 1 -  
Switzerland

Version 11. Revision date 11/23

# extendSURE™

**CONTROL****LOT SPECIFIC VALUE ASSIGNMENT CARD****LEVEL 2**

## Haemoglobin A1c

NGSP ALIGNED

IFCC ALIGNED

**UNITS**

%

mmol/mol

g/dL

mmol/L

**MEAN**

10.1

87

1.29

0.80

**RANGE**

(8.4-11.8)

(73-101)

(1.07-1.51)

(0.66-0.94)

## Total Haemoglobin

g/dL

mmol/L

14.62

9.07

(11.97-17.27)

(7.43-10.71)

**MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY****BECKMAN COULTER AU CHEMISTRY SYSTEMS**

For use with Beckman Coulter HbA1c Reagent Kit Part Number B00389 only

**LOT****4307****CE 2797****EC REP****USE BY 2026-04-30****IVD****CH REP**Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New ZealandEMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The NetherlandsBeckman Coulter International S.A.  
Rue Juste-Olivier 22,  
1260 - Nyon 1 -  
Switzerland

Version 11. Revision date 11/23