

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 1

Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

5.6

(4.6-6.6)

IFCC ALIGNED

mmol/mol

38

(31-45)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC AU CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT

4307

CE 2797



USE BY 2026-04-30



IVD

Ver 04. Revision date 11/23



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EC REP

EMERGO EUROPE
Westervortsedijk 60,
6827 AT Arnhem
The Netherlands

CH REP

Beckman Coulter International S.A.
Rue Juste-Olivier 22,
1260 - Nyon 1 -
Switzerland

extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 2

Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

10.0

(8.3-11.7)

IFCC ALIGNED

mmol/mol

86

(72-100)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER **DxC AU** CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT

4307

CE 2797

EC REP



USE BY 2026-04-30



IVD

CH REP

Ver 04. Revision date 11/23



Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New Zealand

EMERGO EUROPE
Westenvoortsedijk 60,
6827 AT Arnhem
The Netherlands

Beckman Coulter International S.A.
Rue Juste-Olivier 22,
1260 - Nyon 1 -
Switzerland