

# extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 1

## Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

5.7

(4.7 - 6.7)

IFCC ALIGNED

mmol/mol

39

(32 - 46)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

### BECKMAN COULTER **DxC AU** CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009

LOT

4324

CE 2797



USE BY 2027-03-31



IVD

Ver 04. Revision date 11/23



Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New Zealand

EC REP

EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands

CH REP

Beckman Coulter International S.A.  
Rue Juste-Olivier 22,  
1260 - Nyon 1 -  
Switzerland

# extendSURE™

CONTROL

LOT SPECIFIC VALUE ASSIGNMENT CARD

LEVEL 2

## Haemoglobin A1c

UNITS

MEAN

RANGE

NGSP ALIGNED

%

9.9

(8.2 - 11.6)

IFCC ALIGNED

mmol/mol

85

(71-99)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

### BECKMAN COULTER **DxC AU** CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Advanced Reagent Kit Ref. Number B93009



Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New Zealand

LOT

4324

CE 2797

EC REP

EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands



USE BY 2027-03-31



IVD

CH REP

Beckman Coulter International S.A.  
Rue Juste-Olivier 22,  
1260 - Nyon 1 -  
Switzerland

Ver 04. Revision date 11/23